

Ver.24081

OneStep RTPCR Kit

Cat.No. ODR61

Description: The OneStep RTPCR Kit contains all necessary reagents for both reverse transcription and PCR amplification to occur in a single reaction tube. Components for both cDNA synthesis and PCR are combined in a single tube, using gene-specific primers and target RNAs from either total RNA or mRNA. Reverse transcription automatically follows PCR cycling without additional steps. The OneStep RTPCR Kit consists of two major components: Enzyme Mix and 5X Reaction Mix, provided a convenient format for highly sensitive and specific RT-PCR using any RNA template. The Enzyme Mix is for optimal cDNA synthesis and PCR amplification. The enzyme mix uses a mixture of M-MLV Reverse Transcriptase and HotStart *Taq* DNA polymerase in an optimized reaction buffer. The 5X Reaction Mix consists of a proprietary buffer system optimized for reverse transcription and PCR amplification, Mg^{2+} optimized for universal use, deoxyribonucleotide triphosphates, and stabilizers.

Storage

- Store at $-20^{\circ}C$ for 12 months
- Prepare aliquots to avoid multiple freeze-thaw cycles.

Features

Component	Volume
Enzyme Mix	100 μ L
5X Reaction Mix	550 μ L

Recommended PCR mixture and cycling

PCR mixture (Reaction vol. 50 μ L)		Cycle		
5X Reaction Mix	10 μ L	45-55 $^{\circ}C$	15-30 min	$\times 1$
Forward primer (10pmol/ μ L)	1 μ L	94 $^{\circ}C$	5 min	$\times 1$
Reverse primer (10pmol/ μ L)	1 μ L	94 $^{\circ}C$	15 sec	} $\times 35-40$
Template RNA (1pg - 1 μ g)	- μ L	55-60 $^{\circ}C$	30 sec	
Enzyme Mix	2 μ L	68-72 $^{\circ}C$	1 min/kb	
Add D.W to	50 μ L	72 $^{\circ}C$	5-10 min	$\times 1$

Note: You may modify the amount of template, extension time, annealing temperature, and the number of PCR cycle according to the target size, primers T_m and the type of templates for amplification.

Important Parameters

1. RNA

- High quality intact RNA is essential for successful full-length cDNA synthesis.
- RNA should be avoiding of any RNase contamination and aseptic conditions should be maintained.

2. Primers

- Gene specific primers (GSP) are recommended. Use of oligo(dT) or random primers are not recommended as they result in generation of non-specific products in the one-step procedure and the amount of RT-PCR product may be reduced.
- A final primer concentration of 0.2 μ M for each primer is generally optimal. For best results, a primer titration using 0.15-0.5 μ M is recommended. Primers should not be self-complementary or complementary to each other at the 3' ends.

3. Magnesium Concentration

- The magnesium is included in the 5X Reaction Mix at a final concentration of 2.3mM, which works well for most RT-PCR.

4. dNTPs Concentration

- dNTPs are included in the 5X Reaction Mix, which is optimal for most reactions.