

Goat Anti-Rabbit IgG (H+L) (Elab Fluor® 594 conjugated)

OPS1060E

Reactivity	Rabbit	Storage	-20°C, Avoid freeze / thaw cycles
Isotype	IgG	Applications	IF

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Product Information

Immunogen	Rabbit IgG
Purification	Affinity purification
Buffer	PBS with 0.1% Proclin 300, 1% protective protein and 50% glycerol, pH 7.4
Storage	Store at -20°C Valid for 12 months. Avoid freeze/ thaw cycles. Protected from prolonged exposure to light.
Dilution	IF:- 1:100 – 300

Other applications have not been tested. Optimal dilutions/concentrations should be determined by the end user.

Research Use

For research use only, not for use in diagnostic procedures.

Legend

Applications: WB-Western Blot; IHC-Immunohistochemistry; IF-Immunofluorescence; IP-Immunoprecipitation; FC-Flow cytometry; ChIP-Chromatin Immunoprecipitation

Reactivity: H-Human; R-Rat; M-Mouse; Mk-Monkey; Dg-Dog; Ch-Chicken; Hm-Hamster; Rb-Rabbit; Sh- Sheep; Pg-Pig; Z-Zebrafish; X-Xenopus; C-Cow.

Please contact Origin Diagnostics and Research for further assistance

www.originlab.in

info@originlab.in

+91-7736237778