

Ver.23081

1X Phosphate Buffered Saline (PBS)

Cat.No.: PB1020

Description

Phosphate Buffered Saline (PBS) is routinely used as a wash buffer in Western blot and immunoprecipitation (IP) procedures, as well as during sample preparation of cell lysates for use in general molecular biology applications.

Storage

Store at room temperature (15-25°C)

Applications

- Suitable for use in wash buffers for Western blots and IP, as well as cell lysate sample preparation.
- It has been used for the suspension of bacterial cells.
- It has been used during the fixation and immunostaining procedures

Notes

- If a precipitate is formed, gently heat the solution till all the precipitate is dissolved.