



ISO 13485:2016 ISO 9001:2015

Ver.25101

Oligo (dT)₁₈

Cat.No.: OCD106

Description: The oligo (dT)₁₈ primer is a synthetic single-stranded 18-mer oligonucleotide with 5'- and 3'-hydroxyl ends.

Product contents:

Oligo (dT)₁₈ (2.8nmol): Lyophilized

Storage

-20°C

Application

First strand cDNA synthesis.

Preparation

Reconstitute in 56µl nuclease free water to obtain 50µM concentration. Store at -20°C after reconstitution. Make aliquots before storage to avoid repeated freeze-thaw cycles.