

Certificate of Analysis

Product Name : Anti-Human 25 –OH VD Monoclonal Antibody
Catalog Number: 25-OHVD102
Lot Number: 20250409
Immunogen: 25-OH VD-BSA antigen
Form: Liquid
Purity: $\geq 90\%$ (by SDS-PAGE)
Concentration: 2.55 mg/mL
Buffer : 0.01M PBS (pH 7.4)
Preservative: 0.05% Proclin-300
Appearance: Clear or slightly opalescent solution
Application: LF/CLIA/EIA
Storage: Store at $-20 \pm 5^{\circ}\text{C}$ for two years.
Avoid repeated freeze- thaws.

Note: For diagnostic use only. Not intended
For therapeutic use.

Quality Department

Date: 2025.08.04



For further details, contact us at

Origin Diagnostics and Research, XIII/712, Origin Building, Clappana, Karunagappally, Kollam-690525
info@originlab.in | +91 7736957333 | www.originlab.in