

Goat Anti-Mouse IgG -FITC Conjugated

OPS020

Description

Storage	: -20°C, Avoid freeze / thaw cycles
Fluorophore	: FITC
Conjugation	: FITC (E _x /E _m :492nm/520nm)
Buffer	: PBS with 0.02% sodium azide and 50% glycerol, pH 7.4

Applications

Recommended dilutions : 1:20 -1:100 for most applications

Note : Optimal dilutions should be determined by the end user.

Research Use

For research use only, not for use in diagnostic procedures.

Legend

Applications: WB-Western Blot; IHC-Immunohistochemistry; IF-Immunofluorescence; IP-Immunoprecipitation; FC-Flow cytometry; ChIP-Chromatin Immunoprecipitation

Reactivity: H-Human; R-Rat; M-Mouse; Mk-Monkey; Dg-Dog; Ch-Chicken; Hm-Hamster; Rb-Rabbit; Sh- Sheep; Pg-Pig; Z-Zebrafish; X- Xenopus; C-Cow.

Please contact Origin Diagnostics and Research for further assistance

www.originlab.in

info@originlab.in

+91-7736237778